

Instruments (Platform)



What is the Flow Cell

Flow Cell for HiSeq



8 Lane




2 Lane

Flow Cell for MiSeq




Library Preparation: Shotgun



- Adapter ligation
 - T-overhangs
 - Forked structure controls orientation
- Library amplification
 - Few cycles
 - Enrich for correctly-adapted fragments
 - Required to complete adapter structure in some protocols
- Size selection
 - Gel excision, AMPure beads
 - Limit insert size as needed, remove artifacts

What is the adapter?

Figure 1 Sequencing Library after Paired-End Sample Preparation







Figure 1-2-1 Pipeline of paired-end sequencing (www.illumina.com)



What is the SAM file?

- SAM stands for Sequence Alignment/Map format

Col	Field	Type	Regexp/Range	Brief description
1	QNAME	String	[!-?A-~]{1,254}	Query template NAME
2	FLAG	Int	[0,2 ¹⁶ -1]	bitwise FLAG
3	RNAME	String	* [!-()+-<>-~] [!-~]*	Reference sequence NAME
4	POS	Int	[0,2 ³¹ -1]	1-based leftmost mapping POSition
5	MAPQ	Int	[0,2 ⁸ -1]	MAPping Quality
6	CIGAR	String	* ([0-9]+[MIDNSHPX=])+	CIGAR string
7	RNEXT	String	* = [!-()+-<>-~] [!-~]*	Ref. name of the mate/next read
8	PNEXT	Int	[0,2 ³¹ -1]	Position of the mate/next read
9	TLEN	Int	[-2 ³¹ +1,2 ³¹ -1]	observed Template LENgth
10	SEQ	String	* [A-Za-z=.]+	segment SEQUENCE
11	QUAL	String	[!-~]+	ASCII of Phred-scaled base QUALity+33

What is the BAM file

- Binary Alignment/Map

